

## PRIMAFLORA ® VR BIO













## PRIMAFLORA® VR BIO



## **→** INSTRUCTIONS FOR USE

• Rehydrate 500 g of **PRIMAFLORA® VR BIO** in 10 liters of mineral water or non-chlorinated water at 25-30 °C with 0.5 L of MCR.

## --> ADDITIONAL INFORMATION

- Request further technical assistance from your AEB Technical agent.
- Do not keep the solution longer than 45 minutes:
- Double the volume with grape juice to extend the life by 3 hours
- Quadruple the volume with grape juice to extend the life by 12 hours.
- Do not use on sulphitized juices



Store in sealed original package, protected from light, in a dry, odor-free place at a temperature between 4 and 7 °C Do not freeze.

Respect the DDM indicated on the packaging. Use soon after opening

• 500 g bag - 10 kg carton (=  $20 \times 500$  g).

